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» Reduced sloping with the new MLPA buffer «

Main advantages

The new MLPA buffer strongly reduces sloping that is frequently observed when DNA is extracted with a magnetic beads-based system. Sloping is a characteristic pattern with lower peak heights for larger fragments, which can lead to variation or the potential loss of signals. If you used the S4 Sample Stabiliser for this purpose in the past, you may no longer need to add this reagent¹.


Figure: Electropherograms from reactions on the same male DNA extracted with a system based on magnetic beads. Left: Using the old MLPA buffer and Ligase-65. Right: Using the new MLPA buffer and Ligase-65. The new MLPA buffer will not eliminate other causes of sloping.

What will change?

The new MLPA buffer contains the anti-oxidant dithiothreitol (DTT). This counters sloping caused by iron ions in the sample DNA, which is a common impurity in DNA samples extracted with magnetic beads-based systems.

The Ligase-65 enzyme remains unchanged, but the β -mercaptoethanol anti-oxidant in its storage buffer will be replaced by DTT because this is less pungent and less toxic.

The new MLPA buffer and Ligase-65 vials can be recognized by a 'new' sticker that will be placed on the caps until the end of 2019. The MLPA protocol has not changed, but old and new reagents cannot be used together in the same experiment.